



Lunar DPX IQ

The GE Lunar DPX IQ combines the excellent precision and clinical utility of DPX bone densitometers with high resolution imaging. The DPX-IQ provides the speed of entry level fan beam densitometers with better image quality. It employs a reliable pencil beam to achieve the low precision error and low radiation dose needed for the clinical environment.

Product Features

Windows 98 Operating System

- 4.7e Software
- AP Spine
- Femur
- Dual Femur
- Total Body

25536 John R Rd. Madison Heights, MI 48071

Phone: 888-614-9920 Fax: 888-614-9940

www.medplusequipment.com / www.dexaservices.com