

## OEC 6600 MINI C-ARM



The OEC 6600 mini C-arm is a great & portable system for hand and foot surgery applications. The GE OEC 6600 doesn't require big space instead it will meet your space and your budget.

### Specifications:

#### Dimensions

Base Dimensions Height: 163.2cm (64.25 inches)

Width: 69.21cm (27.25 inches)

Depth: 69.21cm (27.25 inches)

#### X-Ray Tube

Type	Fixed anode, tungsten filament
Target Angle	13.5 degrees
Focal Spot	0.25 x 0.25 millimeters
Filtration Inherent	0.2 mm Al equivalent
Total	1.8 mm Al equivalent
Beam Quality Minimum HVL	2.6 mm Al @ 70 kV
Leakage Technique Factors	75 kV and 100 pA
X-ray Tube Rating	80 kilovolts maximum

#### C-Arm Extension

Maximum Front 15.94cm (62.75 inches) Max.

Height Up 17.90cm (705 inches) Mm.

Height Down 22.30cm (8.75 inches)