

# V-MAX

## VASCULAR IMAGING TABLES

*With all the features of the competition...  
At a fraction of the cost*

**4-WAY FLOAT &  
TRENDLENBURG**



**V-MAX**

24" wide top shown

**V-MAX-T**



### STANDARD FEATURES

- Table top automatically stops at level
- Diving board style for easy c-arm access
- Large patient capacity 500lb pound load
- 2" patient high quality comfort pad
- Facial cut out access
- Light weight – easy to physically move the table
- 24" wide top eliminates the needs for arm boards
- 27" wide top is also available
- Easy patient transfer to table
- Low height for patient boarding
- Any OR accessory can be rigidly attached to the table top

**STI**  
SURGICAL TABLES INC.

# Vascular X-Ray tables

## V-MAX

HEIGHT AND 4 WAY FLOAT


- Height 28-38
- Lateral float 14"
- Longitudinal float 36"



## V-MAX-T

HEIGHT, 4 WAY FLOAT AND TRENDELENBURG

- Height 30-40
- Lateral float 14"
- Longitudinal float 36"
- Trendelenburg +/-16 degrees



2" thick high quality patient comfort pad

Carbon fiber top of 24" width eliminates need for armboards and 84" length can accommodate even the tallest patients. 27" wide top also available.

Control Handle for easy table movement

Patient Straps

Facial Access Hole

Bellows to cover and protect table motor components

High end actuator motor can lift over 500 pound load

Easy Lock/Unlock Castors

Diving board style easy c-arm access

Foot Controls

Lightweight which makes it easy for staff to move the table around the room or room to room



For a price quote or for more info call

**888-737-5044**

Or check out our website at [www.surgicaltablesinc.com](http://www.surgicaltablesinc.com)

### OPTIONS ▼


#### JOYSTICK CONTROL



#### MUSHROOM JOYSTICK



#### MODULAR SIDE RAIL SYSTEM



#### SURGICAL ARMBOARDS

(Removeable armboard rotates 180 degrees)



#### PROCEDURAL ARMBOARDS

(Each armboard rotates 180 degrees)



#### LOLLYPOP ARMBOARD



#### RECTANGULAR ARMBOARD



#### VMAX STRAPS can be undone & left hanging

