



EasyOne® Plus

The worldwide leader in spirometry

Easy, Fast, Reliable, Accurate




Features

TrueFlow™ Advantages

- Not influenced by heat, humidity, barometric pressure, or altitude
- No cleaning, no maintenance
- Extremely high accuracy for low flows

Essential Features

- USB cradle for direct connection to printer or PC
- Flexible EMR capabilities
- True portability for testing anytime, anywhere
- Powered by 2 standard AA batteries
- Vast memory
- Selectable predicted values and interpretation
- Automated quality control



System Specifications

Print out: Direct to printer or via PC, plain paper

Patient age: 4 and up

Language: English, French, Spanish, Italian, German

Predicted values: NHANES III, Stanojevic 2009, Knudson 1976, Knudson 1983, Crapo, Morris, Chermiak, ERS (ECCS), and more
 Pediatrics: Polgar, HSU, Dockery, Zapletal

Battery life: Approximately 400 tests

Systems

Part Numbers	Description
1000-1 All-inclusive EasyOne Plus diagnostic spirometry system	Spirometer, USB cradle, printer, software, case of 50 spirettes, 4-year extended warranty
2001-2P EasyOne Plus diagnostic spirometer	Spirometer, USB cradle, printer, software
2001-2NP EasyOne Plus diagnostic spirometer	Spirometer, USB cradle, software

Disposables

Part Numbers	Description
2050-1 spirettes	Case of 50
2050-5 spirettes	Case of 200

Test Modes

- FVC, FVL, SVL, MVV, post testing
- Tidal breathing selectable

Proven Worldwide

- NASA, COPD Gene, BOLD
- Meets ATS/ERS, SSA Disability and NIOSH/OSHA standards.

