

LifeDop™ 250 (Model # L250)



Summit Doppler Systems offers a full line of Doppler systems with the sensitivity and performance you depend on. The **LifeDop 250** large display, hand-held Doppler has excellent sound quality and accepts the full line of non-directional probes for obstetrical or vascular applications. This system is ergonomically designed to fit in the palm of your hand with comfort and allows easy access to each control feature. Each system is individually tested and inspected to ensure the highest quality standards.

Features:

- Large display for indication of heart-rate and battery status
- Superior sound quality from high-fidelity speaker
- Extended battery life
- Enhanced probe sensitivity
- Interchangeable probe capability
- User replaceable cable
- Powered by 3 AA batteries
- Each main unit comes with choice of one probe
- Dimensions (h w l): 140 x 70 x 35 mm
- Weight: 320 grams
- Full 2 year product warranty

DISTRIBUTED BY



SOUTHEASTERN
BIOMEDICAL ASSOCIATES INC

3546 Hickory Blvd. Hudson, NC 28638
888.310.7322 www.sebiomedical.com
Email: info@sebiomedical.com

Summit Doppler Systems, Inc.
4680 Table Mountain Dr. #150
Golden, Colorado 80403
303-423-7572 or 800-554-5090
Fax # 303-940-7165
www.summitdoppler.com