



Wallach Summit Doppler Vantage ABI™ – Simple One-touch ABI Screening in Just 3 Minutes

Higher Throughput ABI Assessment for Your Practice

The Vantage ABI conducts fast, accurate ABI screenings in minutes to assist with the early detection of Peripheral Arterial Disease (P.A.D.). The Vantage ABI completes a single level exam using Digital Fourier Oscillometry (DFO). Summit Doppler worked with researchers from the Cleveland Clinic to create a clinically validated unique DFO Algorithm*.

System features include:

- Speedy throughput with large display and intuitive touchscreen navigation
- Convenient USB drive for data storage, transfer, and EMR interface
- Blood pressure cuffs for obtaining pressures and ankle waveforms

VANABI – Package 1: One-touch, three-minute ABI exams with USB data storage and cuffs included.

VANMAX – Package 2: One-touch, three-minute ABI exams with USB data storage and 4 cuffs, plus built-in printer, basket, and durable stand.



* Patent: WO2012040428A2

Summit Doppler™

a WALLACH brand



| ITEM NUMBER | DESCRIPTION |
|-------------|---------------------------------------|
| VANABI | Vantage ABI – Package 1 |
| VANMAX | Vantage ABI with Stand – Package 2 |
| CUF0002 | 10 cm Blood Pressure Cuff (1) |
| CUF0003 | 12 cm Blood Pressure Cuff (1) |
| K180 | Adhesive Backed Label Paper (5 Rolls) |
| K185 | Label Paper (5 Rolls) |

Summit Doppler Vantage ABI

| | |
|------------------------|---|
| Main Unit Dimensions | 6.4 in. x 11 in. x 7.3 in. (16.25 cm x 27.9 cm x 18.5 cm) |
| Main Unit Weight | 4.6 lbs (2.1 Kg) |
| Mounting | Self-supported (rubber feet) or stand mounted (optional) |
| Portability | Built-in carry handle |
| Cuff Pressure Range | 0 to 280 mmHg |
| Pressure Gage Accuracy | +/-3 mmHg 0–200 mmHg, +/-2% >200 mmHg |
| Deflation Rate | 2.5 mmHg/second nominal |
| Inflation Pressure | 180 mmHg (default), 280 mmHg (user selected) |
| Deflation Pressure | 60 mmHg (default), 35 mmHg (user selected) |
| Pulse Rate Range | 43 to 120 BPM |
| ABI Determination Time | 3 minutes |
| Determination | Pulse Volume Oscillometry |
| Waveform Type | Pneumatic plethysmography (Pulse Volume Recording PVR) |
| Patient Connections | 4 vascular cuffs, 6 ft. hose length (ea.) |
| Clinical Sensitivity | >90% |
| Clinical Specificity | >92% |
| Standards | IEC60601-1, IEC60601-1-2 EMC |
| Classification | Class A |
| Hose/Cable Length | 6 ft. |

Wallach works for you.

Wallach Surgical offers industry-leading medical systems and instruments, including colposcopes, cryosurgical and electrosurgical units, digital and Doppler imaging devices, and a broad range of instruments and disposables that keep practices and hospitals running safely and smoothly.

Physicians have relied on our diverse product portfolio for more than 30 years. Moving forward, we are committed to expanding into new technologies, including our acquisition of Summit Doppler and our ongoing research and development efforts.

WALLACH®

Wallach Surgical Devices
95 Corporate Drive
Trumbull, CT 06611 U.S.A.

1.800.243.2463

www.wallachsurgical.com

© Wallach Surgical Devices 82492 Rev. 4/14