

INDUSTRY LEADING!

UROMAX

X-RAY IMAGING TABLES FROM SURGICAL TABLES INC.

With all
the features of
the competition...
At a fraction of
the cost!

STANDARD FEATURES

- Optional Full Body Extension: Allows table to be easily converted to a full body table for general x-ray imaging
- Remove and install leg extension or full body extension with fluid catcher in place and patient on the table
- Leg extension easy to remove & light weight
- Diving board style for easy c-arm access
- Large patient capacity 500lb pound load
- 2" patient high quality comfort pad
- Light weight – easy to physically move the table
- 24" wide top eliminates the needs for arm boards
- Easy patient transfer to table
- Low height for patient boarding
- Any OR accessory can be rigidly attached to the table top

MEDICAL SPECIALTIES SUPPORTED:

- Pain Management
- Minimally Invasive Spine Surgery
- GI Procedures
- Ultrasound Procedures
- Orthopedic
- Endovascular
- X-Ray Procedures
- Gall Bladder
- Prostate
- Cystoscopy
- Laparoscopy
- Ureteroscopy
- Lithotripsy
- Bladder Access
- Urinary Tract Infections
- Kidney Stone removal
- Percutaneous nephrolithotomy
- Vasectomy



URO top only



OPTIONAL Full body extension
for standard imaging procedures
such as Pain Management



URO Top with
optional foot board

Optional Adjustable Fowler
backrest for Gynecology &
Urodynamics









URO Top with optional full body extension



Tel: 888-737-5044 2 Debusch Ave Unit C3
 Tel: 978-777-4031 Middleton, MA 01949
 Fax: 978-750-0877 USA


www.surgicaltables.com

The URO-MAX Urological series X- Ray tables from Surgical Tables Inc.



Carbon fiber top of 24" width eliminates need for armboards and 84" length (55" URO top and 29" leg extension) can accommodate even the tallest patients.

Optional Stirrups for leg placement



Optional Handset

Bellows to cover and protect table motor components

High end actuator motor can lift over 500 pound load



Foot Controls

Optional Stirrup clamps

Lightweight which makes it easy for staff to move the table around the room or room to room

Optional Fluid Catcher

Diving board style easy c-arm access

Patient Straps

2" thick high quality patient comfort pad

Facial access hole



For a price quote or for more info call

888-737-5044

OPTIONS


LEG STIRRUPS




LEG STIRRUP CLAMP




FULL BODY EXTENSION



FOWLER BACKREST



SURGICAL ARMBOARDS
(Removeable armboard rotates 180°)




LEG EXTENSION - 24" W X 29" L



FLUID CATCHER



PAPER ROLL



www.surgicaltables.com